


*A Breakthrough Design Enhancement!*  
**Paragon CRT® Dual Axis™**




Just watch us®

800-528-8279 | [www.paragoncrt.com](http://www.paragoncrt.com)

**Paragon CRT® Dual Axis™ Will Accommodate Up To 30% More Corneal Reshaping Patients**



Mal-fitting spherical CRT®



Well centered CRT® Dual Axis™

FDA Approved for up to 6D of myopia, with or without 1.75D astigmatism.

- ✘ Paragon CRT® Dual Axis™  
Allows for independent manipulation of a second RZD and LZA.
- ✘ Paragon CRT® Dual Axis™  
Provides a thin-lens design while maintaining flexure resistance.
- ✘ Paragon CRT® Dual Axis™  
Features improved centration through auto-orientation without prism or slab off.
- ✘ Paragon CRT® Dual Axis™  
Is recommended when any of the following are present:
  - Corneal cylinder.
  - 15u elevation difference at 8mm diameter.
  - Incomplete circumferential engagement
  - Decentration of unknown etiology.